


MICS / VATS

*Minimal Invasive Cardiac Surgery
Video Assisted Thoracoscopic Surgery*





Instrument-Tray

Art. RE 1041

for the washing, sterilisation and storage of MICS/VATS instruments


- *stainless steel for washing and careful sterilisation*
- *storage of delicate and fine instruments*


MICS / VATS

Index


Foerster Grasping Clamp	4
Duval & Grasping Clamp	6
Harken & Satinsky Clamp	7
Dissector Sliding Clamp	10
Dissector Sliding Clamp I Debakey & Dennis	11
Sliding Scissors	12
Needle Holder TC/HM Slider	14
Double Action Grasper Clamp Hallmann	15
Double Action Scissors (Metzenbaum)	16
Instrumentation Supplies	
MICS Suction Instrument, Knot Pusher	17
Cardioplegia Needle, Magnet Finder	18

Legend

- 
DEBAKEY
 1x2 rows of microscopic teeth
- 
 2x3 rows of microscopic teeth
- 
Longitudinal teeth
 with longitudinal serrations
- 
Querzahn 60°
- 
Glassman Allis
 5x6 teeth



selected clamps are offered
with / or without ratchet
 to facilitate ease of dissection



All graphics are original size (100%) if not otherwise noted.

Overall length depends on jaw and handle style.

All measurements, pressures and angles are typical and may vary due to production.

Not all products are available in every country.



This Reuchlen catalog is effective with August 2023.

The content is exclusively for information purposes, can be changed without prior notice at any time and can not be seen as our obligation. The August Reuchlen GmbH grants no warranty or guarantee with regard to correctness and accuracy of the content in this catalog.


Instrument design can be subject to change or modification without notice.

© Copyright August Reuchlen GmbH®, Tuttlingen, Germany

FOERSTER GRASPING CLAMP



	RC 1826-07	RC 1830
Art. ratchet shaft	with 7mm	with 10mm
working length	240mm	260mm
total length	335mm	350mm
jaw tip	4mm	5mm



	RC 1859-07	RC 1846-10
Art. ratchet shaft	with 7mm	with 10mm
working length	180mm	250mm
total length	280mm	340mm
jaw tip	5mm	5mm

FOERSTER GRASPING CLAMP




RC 1829



6x12mm
oval ring
external
dimension



Art. ratchet shaft	RC 1827-05 with 5mm	RC 1858-05 with 5mm
working length	240mm	180mm
total length	335mm	280mm
jaw tip	5mm	5mm
	RC 1827-05S	
jaw tip	3mm	



6x12mm
oval ring
external
dimension

7x13mm
oval ring
external
dimension




Art. ratchet shaft	RC 1828-05 with 5mm	RC 1829-07 with 7mm
working length	240mm	240mm
total length	335mm	335mm
jaw tip	4mm	5mm
	RC 1828	RC 1829
shaft	10mm	10mm
jaw tip	5mm	5mm





DUVAL & GRASPING CLAMP




RC 1831



Duval Clamp




				
Art.	RC 1831	RC 1852-10	RC 1832	RC 1850-10
ratchet	with	with	with	with
shaft	10mm	10mm	10mm	10mm
working length	240mm	240mm	245mm	250mm
total length	330mm	330mm	335mm	340mm
jaw width	14mm	14mm	20mm	25mm

Grasping Clamp



		
Art.	RC 1836	RC 1849-10
ratchet	with	with
shaft	10mm	10mm
working length	240mm	240mm
total length	330mm	330mm
jaw length	67mm	67mm
jaw tip	5mm	10mm

HARKEN & SATINSKY CLAMP





Art. ratchet shaft
working length
total length

RC 1796
with
10mm
260mm
349mm

RC 1798
without
10mm
260mm
349mm

RC 1795
with
10mm
260mm
347mm



RC 1797
without
10mm
260mm
347mm




Art. ratchet shaft
working length
total length

RC 1799
with
10mm
270mm
357mm


RC 1838
with
10mm
230mm
316mm




HARKEN & SATINSKY CLAMP



RC 1788




	RC 1786	RC 1787	RC 1788	RC 1789
Art. ratchet shaft	with 10mm	with 10mm	with 10mm	with 10mm
working length	200mm	290mm	285mm	280mm
total length	290mm	380mm	375mm	370mm
jaw length	50mm	115mm	110mm	105mm





Art. ratchet shaft	RC 1811-07 without 7mm
working length	230mm
total length	330mm
jaw length	65mm


HARKEN & SATINSKY CLAMP




RC 1794




	 Art. ratchet shaft working length total length jaw length	RC 1794 with 10mm 200mm 290mm 82mm	RC 1793 with 10mm 240mm 330mm 62mm	RC 1792 with 10mm 260mm 360mm 93mm
--	---	---	---	---

	 Art. ratchet shaft working length total length jaw length	RC 1699 with 10mm 190mm 290mm 82mm	RC 1698 with 10mm 240mm 330mm 62mm	RC 1697 with 10mm 260mm 360mm 92mm
--	---	---	---	---


DISSECTOR SLIDING CLAMP



RC 1833




Art. ratchet shaft	RC 1833 with 10mm	RC 1834 without 10mm	RC 1812-10 without 10mm
working length	240mm	240mm	250mm
total length	330mm	330mm	340mm



Art. ratchet shaft	RC 1813-10 with 10mm	RC 1813-11 without
working length	250mm	
total length	340mm	

DISSECTOR SLIDING CLAMP | DEBAKEY & DENNIS



RC 1835

Dennis Slider




Art. ratchet shaft	RC 1835 with 10mm	RC 1837-07 without 7mm
working length	240mm	230mm
total length	330mm	330mm
jaw length	65mm	65mm
Art. ratchet shaft	RC 1835-07 with 7mm	
working length	230mm	
total length	330mm	
jaw length	65mm	



DeBakey Slider



Art. ratchet shaft	RC 1851-05 without 5mm
working length	190mm
total length	290mm
Art. ratchet	RC 1851-06 with



SLIDING SCISSORS



RA 8110




Art.	RA 8110	RA 8112
ratchet	without	without
shaft	10mm	10mm
working length	250mm	250mm
total length	360mm	360mm
jaw length	55mm	55mm



Art.	RA 8114-07	RA 8114-360	RA 8114
ratchet	without	without	without
shaft	7mm	10mm	10mm
working length	165mm	250mm	190mm
total length	290mm	360mm	290mm
jaw length	35mm	35mm	35mm

SLIDING SCISSORS




RA 8116



	RA 8116	RA 8118-05	RA 8120-05
Art. ratchet	without	without	without
shaft	5mm	5mm	5mm
working length	228mm	195mm	165mm
total length	355mm	320mm	290mm
jaw length	35mm	35mm	30mm


NEEDLE HOLDER TC/HM SLIDER



RN 2732

TC

Art.
ratchet
shaft
working length
total length
suture size max.
jaw tip width



RN 2730
with
10mm
200mm
285mm
4-0 thread
3mm

RN 2732
with
10mm
200mm
285mm
4-0 thread
3mm

TC

Art.
ratchet
shaft
working length
total length
suture size max.
jaw tip width



RN 2740
with
10mm
180mm
285mm
2-0 thread
4mm


RN 2741-07
with
7mm
240mm
340mm
2-0 thread
4mm

shaft

RN 2741
10mm

TC

Art.
ratchet
shaft
working length
total length
suture size max.
jaw tip width




RN 2742
with
10mm
180mm
285mm
2-0 thread
4mm

RN 2743-07
with
7mm
240mm
340mm
2-0 thread
4mm

DOUBLE ACTION GRASPER CLAMP HALLMANN




RC 4493




	RC 4493	RC 4494
Art. ratchet shaft	with 8mm	with 8mm
working length	151mm	168mm
total length	345mm	400mm
jaw	with ring tip cross teeth	with ring tip cross teeth

DOUBLE ACTION SCISSORS




RA 8104



	RA 8100	RA 8104
Art. shaft	10mm	10mm
working length	214mm	214mm
total length	326mm	366mm
jaw length	73mm	73mm



RA 8102



	RA 8102	RA 8106
Art. shaft	10mm	10mm
working length	208mm	208mm
total length	320mm	360mm
jaw length	67mm	67mm

SUPPLIES



Art. **RL 2021**
shaft 5mm
working length 350mm
total length 475mm
jaw tip width 5mm

MICS Suction Instrument
RL 2021



Art. **RC 0451**
shaft 5mm
working length 230mm
total length 330mm
jaw tip length 9mm

Knot Pusher
RC 0451

SUPPLIES



Art. **RC 0462**
shaft 5mm
working length 310mm
total length 340mm

Cardioplegia Needle
RC 0462



Art. **RC 0457**
shaft 5mm
working length 240mm
total length 340mm
jaw length 20mm,
Ø10mm

Magnet Finder
RC 0457





REU-D / MV-V2 / 2023-08

80 YEARS
1943-2023



August Reuchlen GmbH

Junkersstr. 8
78532 Tuttlingen
Germany
Tel.: +49 (0) 74 61 / 96 52 - 0
Fax: +49 (0) 74 61 / 96 52 - 50
info@reuchlen.com

www.reuchlen.com

